

SPECIFICATIONS

Main Function	<ol style="list-style-type: none"> 1. A professional salon spa-treatment at home 10 minutes a day promotes youthful skin. 2. Facial massager combination of negative ionizer, positive ionizer, 5MHz ultrasonic all in one beauty device. 3. Enhance the look of skin affected by conditions like poor skin texture and uneven pigmentation. 4. Ultrasonic and negative ions treats your skin for better absorption of skin care products. 5. Positive ions cleanse skin by causing dead skin cells, dirt, excess oil, etc., to cling to each other and form particles for easier removal from skin surface.
Advantage	<ol style="list-style-type: none"> 1. Multi-function Ultrasonic, negative ions and positive ions combined in one machine; 2. High Ultrasonic Frequency at 5MHz specially designed for face usage which is more safe than other in 1MHz or 3MHz; 3. Cordless with rechargeable battery built-in, more convenient for use at anywhere; 4. USB charger convenient for charging during traveling; 5. High end elegant industrial design with pearl white painting; 6. 10 minutes automatically power off.
Accessories	USB line Wall adaptor
Package	Gift box size: 182x135x54mm Master Carton: @36pcs/6.8kgs/8.0kgs/422*342*386mm
Warranty	One (1) Year from the date of purchase from any defect caused by material or workmanship

Includes:

- The Bella Leeza™ Ultrasonic +ION Face Massager
- AC Adapter

Bella Leeza
P.O. Box 76
Glenview, IL 60025-0076
USA
tel: +1 847-465-1100

The Bella Leeza™ Ultrasonic +ION Face Massager™

OPERATION MANUAL

(Please read the manual carefully before use)

How Does It Work?

The sonic vibrations and ION technology action generated by **The Bella Leeza™ Ultrasonic+ION Face Massager** helps revive cell tissue, excrete skin wastes, and facilitate skin metabolism.

The application of cosmetic ultrasound is commonly used by spa and salon professionals to “plump up” tired, aging skin and smooth away the appearance of fine lines and wrinkles. Ultrasonic treatment also can also help shrink enlarged pores, relieve puffy eyes, reduce the appearance of cellulite, fade dark eye circles, freckles and old age spots, improve the appearance of acne and even repair past scar damage. The process is safe, painless, gradual and gentle and appropriate for all skin types.

The Bella Leeza™ Ultrasonic +ION Face Massager uses high level sound wave technology to penetrate below the surface of the skin promoting cellular renewal and repair, muscle toning, increased blood circulation, lymphatic drainage and a reduction in puffiness and swelling. Super ION releasing technology, promotes increased blood circulation, accelerates metabolism, makes skin cells regain vitality, resulting in decreased number and depth of wrinkles, as well as prevents new wrinkles from forming.

Maximizes the Effectiveness of Your Existing Skin Care Products

The overall effectiveness of your favorite creams, serums and lotions will be vastly improved with treatment. The ION and Sonic waves emitted by **The Bella Leeza™ Ultrasonic +ION Face Massager** make the skin more receptive and allow applied topical to penetrate much deeper below the surface of the skin.

Results Are Gradual And Can Vary By Individual Skin Type

Those using ultrasonic facial treatment for anti-aging purposes may not realize the full benefits of treatment for several weeks as treatment will be gradual. Ultrasonic treatment has been shown to produce an immediate and temporary lifting effect, while continued daily application will provide cumulative, more long lasting results.



Benefits Of Ultrasonic Facial Treatment

- Softening of Fine Lines
- Relief of Mild Acne
- Refinement of Large Pores
- Firming and Toning of Sagging Facial Skin
- Decongestion of Puffy Eyes
- Fading of Dark Eye Circles
- Reduction in The Appearance of Cellulite
- Renewal of Skin Texture and Tone
- Fading of Age Spots and Freckles Due To Sun Damage
- Improvement In The Appearance of Scars and Stretch Marks

Features

- Easy-To-Use And Ideal For Home Use.
- Large Treatment Head For Improved Effectiveness of Sonic and ION Technology
- Offers The Benefit of Portability For Frequent Travelers.
- A Simple and Effective Treatment Requiring Only Minutes Per Day.
- Effective And Safe For All Skin Types And Ethnicities.

General Precautions

- Do Not Use If Pregnant, Have a Pacemaker or Medical Condition.
- Always Keep Silver Massage Head in Contact With Skin.
- Never Use Internally (i.e. in Mouth, Nose, etc.)
- Avoid Using Near The Mouth, Ears and Center Throat Area.

How often and how long should I use the System?

The Ultrasonic Home Facial Massager may be used daily. Total treatment time should not exceed 5 minutes per day for the facial area. Sensitive skin types should limit treatment time to 2-3 minutes per day. Treatment of general body area including arms, legs, buttocks, thighs and stomach should not exceed 10 minutes per day.

How long will it take to see results?

Everyone's skin varies, but most people feel and see a difference in about 3 weeks. Because of differences in skin conditions, some people will see the benefits almost immediately and others will need more treatments. Typically within a few uses, people feel a tingling, slight warmth and a tightening of their skin. Because treatments are progressive, we recommend treatments for at least ten (10) weeks, for best results.

Is home ultrasonic facial treatment safe?

Ultrasonic facial treatment is considered a safe, easy and gentle skin rejuvenation technology however the following precautions should always be taken:

- Do Not Use If Pregnant, Have a Pacemaker or Medical Condition.
- Always Keep Silver Massage Head in Contact With Skin.
- Never Use Internally (i.e. in Mouth, Nose, etc.)
- Avoid Using Near The Mouth, Ears and Center Throat Area.

Why can't I hear The Bella Leeza™ Ultrasonic +ION Face Massager when switched on?

The sound waves resonate well above our hearing level and move too quickly for our nerve endings to register. As such, this high frequency vibration is often not readily visible or felt. You can test the invisible ultrasonic wave action by placing a small drop of water directly on the wand's silver massage head. When the device is switched on, you should immediately notice vibration-type patterns on the drop of water. This means that your Ultrasonic System is indeed working!

How do you use the system?

1. Charge the device for at least three hours prior to using it.
2. Thoroughly cleanse and dry your skin prior to treatment. The elimination of dead skin cells is the first step in any effective facial treatment program.
3. If using skincare products, apply generously immediately prior to treatment.
4. Turn the unit on and choose the mode depending on the goal of the treatment.
5. Start on **LOW** and after few treatments if no irritation occurs, start using **HIGH**.
6. The unit will shut off automatically after 10 minutes.

OPERATION



1. **General treatment:** around 4 minutes for each side of face, overall 8 minutes.



2. Follow application directions on both sides of your face. At each application direction, repeat the movement 8 times before going to the next direction.



3. For enhanced results, please continue to use your skin care products.